

7 Literaturverzeichnis

1. Handwerker, H. O. *Einführung in die Pathophysiologie des Schmerzes* (Springer Verlag, Berlin, Heidelberg, 1999).
2. Kandel, E. R., Schwartz, J. H. & Jessel, T. M. *Principles of neural science* (1991).
3. Julius, D. & Basbaum, A. I. Molecular mechanisms of nociception. *Nature* **413**, 203-10 (2001).
4. Fields, H. *Pain* (McGraw-Hill, New York, 1987).
5. Scholz, J. & Woolf, C. J. Can we conquer pain? *Nat Neurosci* **5 Suppl**, 1062-7 (2002).
6. Hunt, S. P. & Rossi, J. Peptide- and non-peptide-containing unmyelinated primary afferents: the parallel processing of nociceptive information. *Philos Trans R Soc Lond B Biol Sci* **308**, 283-9 (1985).
7. Snider, W. D. & McMahon, S. B. Tackling pain at the source: new ideas about nociceptors. *Neuron* **20**, 629-32 (1998).
8. Sah, D. W., Ossipo, M. H. & Porreca, F. Neurotrophic factors as novel therapeutics for neuropathic pain. *Nat Rev Drug Discov* **2**, 460-72 (2003).
9. Cesare, P. & McNaughton, P. Peripheral pain mechanisms. *Curr Opin Neurobiol* **7**, 493-9 (1997).
10. Woolf, C. J. & Salter, M. W. Neuronal plasticity: increasing the gain in pain. *Science* **288**, 1765-9 (2000).
11. Melzack, R. & Wall, P. D. Pain mechanisms: a new theory. *Science* **150**, 971-9 (1965).
12. Bourne, H. R. How receptors talk to trimeric G proteins. *Curr Opin Cell Biol* **9**, 134-42 (1997).
13. Gudermann, T., Kalkbrenner, F. & Schultz, G. Diversity and selectivity of receptor-G protein interaction. *Annu Rev Pharmacol Toxicol* **36**, 429-59 (1996).
14. Jordt, S. E., Tominaga, M. & Julius, D. Acid potentiation of the capsaicin receptor determined by a key extracellular site. *Proc Natl Acad Sci U S A* **97**, 8134-9 (2000).
15. Steen, K. H. & Reeh, P. W. Actions of cholinergic agonists and antagonists on sensory nerve endings in rat skin, in vitro. *J Neurophysiol* **70**, 397-405 (1993).
16. Moriyama, T. et al. Possible involvement of P2Y2 metabotropic receptors in ATP-induced transient receptor potential vanilloid receptor 1-mediated thermal hypersensitivity. *J Neurosci* **23**, 6058-62 (2003).
17. Pin, J. P. & Bockaert, J. Get receptive to metabotropic glutamate receptors. *Curr Opin Neurobiol* **5**, 342-9 (1995).
18. Chuang, H. H. et al. Bradykinin and nerve growth factor release the capsaicin receptor from PtdIns(4,5)P2-mediated inhibition. *Nature* **411**, 957-62 (2001).
19. Conner, A. C. et al. Interaction of calcitonin-gene-related peptide with its receptors. *Biochem Soc Trans* **30**, 451-5 (2002).
20. Mullins, D. E., Zhang, X. & Hawes, B. E. Activation of extracellular signal regulated protein kinase by neuropeptide Y and pancreatic polypeptide in

- CHO cells expressing the NPY Y(1), Y(2), Y(4) and Y(5) receptor subtypes. *Regul Pept* **105**, 65-73 (2002).
21. Przewlocki, R. & Przewlocka, B. Opioids in chronic pain. *Eur J Pharmacol* **429**, 79-91 (2001).
 22. Vanegas, H. & Schaible, H. G. Prostaglandins and cyclooxygenases [correction of cyclooxygenases] in the spinal cord. *Prog Neurobiol* **64**, 327-63 (2001).
 23. Di Marzo, V., De Petrocellis, L., Fezza, F., Ligresti, A. & Bisogno, T. Anandamide receptors. *Prostaglandins Leukot Essent Fatty Acids* **66**, 377-91 (2002).
 24. Clapham, D. E. Calcium signaling. *Cell* **80**, 259-68 (1995).
 25. Berridge, M. J. Inositol trisphosphate and calcium signalling. *Nature* **361**, 315-25 (1993).
 26. Tsien, R. W., Lipscombe, D., Madison, D., Bley, K. & Fox, A. Reflections on Ca(2+)-channel diversity, 1988-1994. *Trends Neurosci* **18**, 52-4 (1995).
 27. Jones, S. & Yakel, J. L. Ca²⁺ influx through voltage-gated Ca²⁺ channels regulates 5-HT₃ receptor channel desensitization in rat glioma x mouse neuroblastoma hybrid NG108-15 cells. *J Physiol* **510** (Pt 2), 361-70 (1998).
 28. Kato, A., Ohkubo, T. & Kitamura, K. Algogen-specific pain processing in mouse spinal cord: differential involvement of voltage-dependent Ca(2+) channels in synaptic transmission. *Br J Pharmacol* **135**, 1336-42 (2002).
 29. Nebe, J., Vanegas, H. & Schaible, H. G. Spinal application of omega-conotoxin GVIA, an N-type calcium channel antagonist, attenuates enhancement of dorsal spinal neuronal responses caused by intra-articular injection of mustard oil in the rat. *Exp Brain Res* **120**, 61-9 (1998).
 30. Bowersox, S. S. et al. Selective N-type neuronal voltage-sensitive calcium channel blocker, SNX-111, produces spinal antinociception in rat models of acute, persistent and neuropathic pain. *J Pharmacol Exp Ther* **279**, 1243-9 (1996).
 31. Atanassoff, P. G. et al. Ziconotide, a new N-type calcium channel blocker, administered intrathecally for acute postoperative pain. *Reg Anesth Pain Med* **25**, 274-8 (2000).
 32. Wermeling, D. et al. Pharmacokinetics and pharmacodynamics of intrathecal ziconotide in chronic pain patients. *J Clin Pharmacol* **43**, 624-36 (2003).
 33. Staats, P. S. et al. Intrathecal ziconotide in the treatment of refractory pain in patients with cancer or AIDS: a randomized controlled trial. *Jama* **291**, 63-70 (2004).
 34. Gullapalli, S. & Ramarao, P. L-type Ca²⁺ channel modulation by dihydropyridines potentiates kappa-opioid receptor agonist induced acute analgesia and inhibits development of tolerance in rats. *Neuropharmacology* **42**, 467-75 (2002).
 35. Barnard, E. A. Receptor classes and the transmitter-gated ion channels. *Trends Biochem Sci* **17**, 368-74 (1992).
 36. Unwin, N. Neurotransmitter action: opening of ligand-gated ion channels. *Cell* **72 Suppl**, 31-41 (1993).
 37. Neher, E. & Sakmann, B. Single-channel currents recorded from membrane of denervated frog muscle fibres. *Nature* **260**, 799-802 (1976).

38. Casteels, R., Raeymaekers, L., Suzuki, H. & Van Eldere, J. Tension response and ⁴⁵Ca release in vascular smooth muscle incubated in Ca-free solution. *Pflugers Arch* **392**, 139-45 (1981).
39. Putney, J. W., Jr. A model for receptor-regulated calcium entry. *Cell Calcium* **7**, 1-12 (1986).
40. Putney, J. W., Jr. & Bird, G. S. The inositol phosphate-calcium signaling system in nonexcitable cells. *Endocr Rev* **14**, 610-31 (1993).
41. Tsiokas, L. et al. Specific association of the gene product of PKD2 with the TRPC1 channel. *Proc Natl Acad Sci U S A* **96**, 3934-9 (1999).
42. Hardie, R. C. TRP channels in Drosophila photoreceptors: the lipid connection. *Cell Calcium* **33**, 385-93 (2003).
43. Minke, B. & Cook, B. TRP channel proteins and signal transduction. *Physiol Rev* **82**, 429-72 (2002).
44. Voets, T. & Nilius, B. The pore of TRP channels: trivial or neglected? *Cell Calcium* **33**, 299-302 (2003).
45. Hotta, Y. & Benzer, S. Abnormal electroretinograms in visual mutants of Drosophila. *Nature* **222**, 354-6 (1969).
46. Montell, C. & Rubin, G. M. Molecular characterization of the Drosophila trp locus: a putative integral membrane protein required for phototransduction. *Neuron* **2**, 1313-23 (1989).
47. Harteneck, C., Plant, T. D. & Schultz, G. From worm to man: three subfamilies of TRP channels. *Trends Neurosci* **23**, 159-66 (2000).
48. Montell, C. et al. A unified nomenclature for the superfamily of TRP cation channels. *Mol Cell* **9**, 229-31 (2002).
49. Clapham, D. E. TRP channels as cellular sensors. *Nature* **426**, 517-24 (2003).
50. Chen, X. Z. et al. Polycystin-L is a calcium-regulated cation channel permeable to calcium ions. *Nature* **401**, 383-6 (1999).
51. Yellen, G. The voltage-gated potassium channels and their relatives. *Nature* **419**, 35-42 (2002).
52. Montell, C., Birnbaumer, L. & Flockerzi, V. The TRP channels, a remarkably functional family. *Cell* **108**, 595-8 (2002).
53. Patapoutian, A., Peier, A. M., Story, G. M. & Viswanath, V. ThermoTRP channels and beyond: mechanisms of temperature sensation. *Nat Rev Neurosci* **4**, 529-39 (2003).
54. Holzer, P. Capsaicin as a tool for studying sensory neuron functions. *Adv Exp Med Biol* **298**, 3-16 (1991).
55. Jancso, G. & Wollemann, M. The effect of capsaicin on the adenylate cyclase activity of rat brain. *Brain Res* **123**, 323-9 (1977).
56. Jancso, N., Jancso-Gabor, A. & Szolcsanyi, J. Direct evidence for neurogenic inflammation and its prevention by denervation and by pretreatment with capsaicin. *Br J Pharmacol* **31**, 138-51 (1967).
57. Caterina, M. J. et al. The capsaicin receptor: a heat-activated ion channel in the pain pathway. *Nature* **389**, 816-24 (1997).
58. Tominaga, M. & Caterina, M. J. Thermosensation and pain. *J Neurobiol* **61**, 3-12 (2004).
59. Hensel, H. & Iggo, A. Analysis of cutaneous warm and cold fibres in primates. *Pflugers Arch* **329**, 1-8 (1971).
60. Iggo, A. Cutaneous thermoreceptors in primates and sub-primates. *J Physiol* **200**, 403-30 (1969).

61. Hensel, H. [Cutaneous heat receptors in primates]. *Pflugers Arch* **313**, 150-2 (1969).
62. Welch, J. M., Simon, S. A. & Reinhart, P. H. The activation mechanism of rat vanilloid receptor 1 by capsaicin involves the pore domain and differs from the activation by either acid or heat. *Proc Natl Acad Sci U S A* **97**, 13889-94 (2000).
63. Tominaga, M. et al. The cloned capsaicin receptor integrates multiple pain-producing stimuli. *Neuron* **21**, 531-43 (1998).
64. Jordt, S. E. & Julius, D. Molecular basis for species-specific sensitivity to "hot" chili peppers. *Cell* **108**, 421-30 (2002).
65. Davis, J. B. et al. Vanilloid receptor-1 is essential for inflammatory thermal hyperalgesia. *Nature* **405**, 183-7 (2000).
66. Caterina, M. J. et al. Impaired nociception and pain sensation in mice lacking the capsaicin receptor. *Science* **288**, 306-13 (2000).
67. Tillman, D. B., Treede, R. D., Meyer, R. A. & Campbell, J. N. Response of C fibre nociceptors in the anaesthetized monkey to heat stimuli: estimates of receptor depth and threshold. *J Physiol* **485 (Pt 3)**, 753-65 (1995).
68. Hu, H. J., Bhawe, G. & Gereau, R. W. t. Prostaglandin and protein kinase A-dependent modulation of vanilloid receptor function by metabotropic glutamate receptor 5: potential mechanism for thermal hyperalgesia. *J Neurosci* **22**, 7444-52 (2002).
69. Bhawe, G. et al. cAMP-dependent protein kinase regulates desensitization of the capsaicin receptor (VR1) by direct phosphorylation. *Neuron* **35**, 721-31 (2002).
70. Hwang, S. W. et al. Direct activation of capsaicin receptors by products of lipoxygenases: endogenous capsaicin-like substances. *Proc Natl Acad Sci U S A* **97**, 6155-60 (2000).
71. Tominaga, M., Wada, M. & Masu, M. Potentiation of capsaicin receptor activity by metabotropic ATP receptors as a possible mechanism for ATP-evoked pain and hyperalgesia. *Proc Natl Acad Sci U S A* **98**, 6951-6 (2001).
72. Zygmunt, P. M. et al. Vanilloid receptors on sensory nerves mediate the vasodilator action of anandamide. *Nature* **400**, 452-7 (1999).
73. Premkumar, L. S. & Ahern, G. P. Induction of vanilloid receptor channel activity by protein kinase C. *Nature* **408**, 985-90 (2000).
74. Winston, J., Toma, H., Shenoy, M. & Pasricha, P. J. Nerve growth factor regulates VR-1 mRNA levels in cultures of adult dorsal root ganglion neurons. *Pain* **89**, 181-6 (2001).
75. Ji, R. R., Samad, T. A., Jin, S. X., Schmall, R. & Woolf, C. J. p38 MAPK activation by NGF in primary sensory neurons after inflammation increases TRPV1 levels and maintains heat hyperalgesia. *Neuron* **36**, 57-68 (2002).
76. Gunthorpe, M. J., Benham, C. D., Randall, A. & Davis, J. B. The diversity in the vanilloid (TRPV) receptor family of ion channels. *Trends Pharmacol Sci* **23**, 183-91 (2002).
77. Prescott, E. D. & Julius, D. A modular PIP2 binding site as a determinant of capsaicin receptor sensitivity. *Science* **300**, 1284-8 (2003).
78. Ahluwalia, J., Yaqoob, M., Urban, L., Bevan, S. & Nagy, I. Activation of capsaicin-sensitive primary sensory neurones induces anandamide production and release. *J Neurochem* **84**, 585-91 (2003).

79. Di Marzo, V. et al. Interactions between synthetic vanilloids and the endogenous cannabinoid system. *FEBS Lett* **436**, 449-54 (1998).
80. Di Marzo, V., Blumberg, P. M. & Szallasi, A. Endovanilloid signaling in pain. *Curr Opin Neurobiol* **12**, 372-9 (2002).
81. Huang, S. M. et al. An endogenous capsaicin-like substance with high potency at recombinant and native vanilloid VR1 receptors. *Proc Natl Acad Sci U S A* **99**, 8400-5 (2002).
82. Chu, C. J. et al. N-oleoyldopamine: A novel endogenous capsaicin-like lipid that produces hyperalgesia. *J Biol Chem* (2003).
83. Sanchez, J. F., Krause, J. E. & Cortright, D. N. The distribution and regulation of vanilloid receptor VR1 and VR1 5' splice variant RNA expression in rat. *Neuroscience* **107**, 373-81 (2001).
84. Mezey, E. et al. Distribution of mRNA for vanilloid receptor subtype 1 (VR1), and VR1-like immunoreactivity, in the central nervous system of the rat and human. *Proc Natl Acad Sci U S A* **97**, 3655-60 (2000).
85. Yiangou, Y. et al. Capsaicin receptor VR1 and ATP-gated ion channel P2X3 in human urinary bladder. *BJU Int* **87**, 774-9 (2001).
86. Van der Aa, F., Roskams, T., Blyweert, W. & De Ridder, D. Interstitial cells in the human prostate: a new therapeutic target? *Prostate* **56**, 250-5 (2003).
87. Kollarik, M. & Udem, B. J. Activation of bronchopulmonary vagal afferent nerves with bradykinin, acid and vanilloid receptor agonists in wild-type and TRPV1-/- mice. *J Physiol* **555**, 115-23 (2004).
88. Ward, S. M., Bayguinov, J., Won, K. J., Grundy, D. & Berthoud, H. R. Distribution of the vanilloid receptor (VR1) in the gastrointestinal tract. *J Comp Neurol* **465**, 121-35 (2003).
89. Renton, T., Yiangou, Y., Baecker, P. A., Ford, A. P. & Anand, P. Capsaicin receptor VR1 and ATP purinoceptor P2X3 in painful and nonpainful human tooth pulp. *J Orofac Pain* **17**, 245-50 (2003).
90. Southall, M. D. et al. Activation of epidermal vanilloid receptor-1 induces release of proinflammatory mediators in human keratinocytes. *J Pharmacol Exp Ther* **304**, 217-22 (2003).
91. Birder, L. A. et al. Vanilloid receptor expression suggests a sensory role for urinary bladder epithelial cells. *Proc Natl Acad Sci U S A* **98**, 13396-401 (2001).
92. Reilly, C. A. et al. Capsaicinoids Cause Inflammation and Epithelial Cell Death through Activation of Vanilloid Receptors. *Toxicol Sci* **73**, 170-81 (2003).
93. Heiner, I. et al. Expression profile of the transient receptor potential (TRP) family in neutrophil granulocytes: evidence for currents through LTRPC2 induced by ADP-ribose and NAD. *Biochem J Pt* (2003).
94. Chen, C. W. et al. Signal transduction for inhibition of inducible nitric oxide synthase and cyclooxygenase-2 induction by capsaicin and related analogs in macrophages. *Br J Pharmacol* **140**, 1077-87 (2003).
95. Birder, L. A. et al. Altered urinary bladder function in mice lacking the vanilloid receptor TRPV1. *Nat Neurosci* **5**, 856-60 (2002).
96. Caterina, M. J., Rosen, T. A., Tominaga, M., Brake, A. J. & Julius, D. A capsaicin-receptor homologue with a high threshold for noxious heat. *Nature* **398**, 436-41 (1999).

97. Kanzaki, M. et al. Translocation of a calcium-permeable cation channel induced by insulin-like growth factor-I. *Nat Cell Biol* **1**, 165-70 (1999).
98. Lewinter, R. D., Skinner, K., Julius, D. & Basbaum, A. I. Immunoreactive TRPV-2 (VRL-1), a capsaicin receptor homolog, in the spinal cord of the rat. *J Comp Neurol* **470**, 400-8 (2004).
99. Ichikawa, H. & Sugimoto, T. Vanilloid receptor 1-like receptor-immunoreactive primary sensory neurons in the rat trigeminal nervous system. *Neuroscience* **101**, 719-25 (2000).
100. Stenholm, E., Bongenhielm, U., Ahlquist, M. & Fried, K. VRL- and VRL-I-like immunoreactivity in normal and injured trigeminal dental primary sensory neurons of the rat. *Acta Odontol Scand* **60**, 72-9 (2002).
101. Bezzerides, V. J., Ramsey, I. S., Kotecha, S., Greka, A. & Clapham, D. E. Rapid vesicular translocation and insertion of TRP channels. *Nat Cell Biol* **6**, 709-20 (2004).
102. Xu, H. et al. TRPV3 is a calcium-permeable temperature-sensitive cation channel. *Nature* **418**, 181-6 (2002).
103. Smith, G. D. et al. TRPV3 is a temperature-sensitive vanilloid receptor-like protein. *Nature* **418**, 186-90 (2002).
104. Peier, A. M. et al. A heat-sensitive TRP channel expressed in keratinocytes. *Science* **296**, 2046-9 (2002).
105. Wissenbach, U., Boddling, M., Freichel, M. & Flockerzi, V. Trp12, a novel Trp related protein from kidney. *FEBS Lett* **485**, 127-34 (2000).
106. Liedtke, W. et al. Vanilloid receptor-related osmotically activated channel (VR-OAC), a candidate vertebrate osmoreceptor. *Cell* **103**, 525-35 (2000).
107. Guler, A. D. et al. Heat-evoked activation of the ion channel, TRPV4. *J Neurosci* **22**, 6408-14 (2002).
108. Watanabe, H. et al. Heat-evoked activation of TRPV4 channels in a HEK293 cell expression system and in native mouse aorta endothelial cells. *J Biol Chem* **277**, 47044-51 (2002).
109. Hensel, H. Thermoreception and temperature regulation. *Monogr Physiol Soc* **38**, 1-321 (1981).
110. Spray, D. C. Cutaneous temperature receptors. *Annu Rev Physiol* **48**, 625-38 (1986).
111. Campero, M., Serra, J. & Ochoa, J. L. C-polymodal nociceptors activated by noxious low temperature in human skin. *J Physiol* **497** (Pt 2), 565-72 (1996).
112. Peier, A. M. et al. A TRP channel that senses cold stimuli and menthol. *Cell* **108**, 705-15 (2002).
113. McKemy, D. D., Neuhausser, W. M. & Julius, D. Identification of a cold receptor reveals a general role for TRP channels in thermosensation. *Nature* **416**, 52-8 (2002).
114. Nealen, M. L., Gold, M. S., Thut, P. D. & Caterina, M. J. TRPM8 mRNA is expressed in a subset of cold-responsive trigeminal neurons from rat. *J Neurophysiol* **90**, 515-20 (2003).
115. Tsavaler, L., Shapero, M. H., Morkowski, S. & Laus, R. Trp-p8, a novel prostate-specific gene, is up-regulated in prostate cancer and other malignancies and shares high homology with transient receptor potential calcium channel proteins. *Cancer Res* **61**, 3760-9 (2001).
116. Schenker, T. & Trueb, B. Down-regulated proteins of mesenchymal tumor cells. *Exp Cell Res* **239**, 161-8 (1998).

117. Jaquemar, D., Schenker, T. & Trueb, B. An ankyrin-like protein with transmembrane domains is specifically lost after oncogenic transformation of human fibroblasts. *J Biol Chem* **274**, 7325-33 (1999).
118. Jordt, S. E. et al. Mustard oils and cannabinoids excite sensory nerve fibres through the TRP channel ANKTM1. *Nature* **427**, 260-5 (2004).
119. Bandell, M. et al. Noxious cold ion channel TRPA1 is activated by pungent compounds and bradykinin. *Neuron* **41**, 849-57 (2004).
120. Kennedy, M. B. Signal-processing machines at the postsynaptic density. *Science* **290**, 750-4 (2000).
121. Kennedy, M. B. The postsynaptic density at glutamatergic synapses. *Trends Neurosci* **20**, 264-8 (1997).
122. Gray, E. G. & Whittaker, V. P. The isolation of nerve endings from brain: an electron-microscopic study of cell fragments derived by homogenization and centrifugation. *J Anat* **96**, 79-88 (1962).
123. tom Dieck, S. et al. Bassoon, a novel zinc-finger CAG/glutamine-repeat protein selectively localized at the active zone of presynaptic nerve terminals. *J Cell Biol* **142**, 499-509 (1998).
124. Hartinger, J., Stenius, K., Hogemann, D. & Jahn, R. 16-BAC/SDS-PAGE: a two-dimensional gel electrophoresis system suitable for the separation of integral membrane proteins. *Anal Biochem* **240**, 126-33 (1996).
125. Dreger, M., Bengtsson, L., Schoneberg, T., Otto, H. & Hucho, F. Nuclear envelope proteomics: novel integral membrane proteins of the inner nuclear membrane. *Proc Natl Acad Sci U S A* **98**, 11943-8 (2001).
126. Klose, J. Protein mapping by combined isoelectric focusing and electrophoresis of mouse tissues. A novel approach to testing for induced point mutations in mammals. *Humangenetik* **26**, 231-43 (1975).
127. O'Farrell, P. H. High resolution two-dimensional electrophoresis of proteins. *J Biol Chem* **250**, 4007-21 (1975).
128. Tsutsumi, S. et al. Propofol activates vanilloid receptor channels expressed in human embryonic kidney 293 cells. *Neurosci Lett* **312**, 45-9 (2001).
129. Schumacher, M. A., Moff, I., Sudanagunta, S. P. & Levine, J. D. Molecular cloning of an N-terminal splice variant of the capsaicin receptor. Loss of N-terminal domain suggests functional divergence among capsaicin receptor subtypes. *J Biol Chem* **275**, 2756-62 (2000).
130. Suzuki, M., Sato, J., Kutsuwada, K., Ooki, G. & Imai, M. Cloning of a stretch-inhibitable nonselective cation channel. *J Biol Chem* **274**, 6330-5 (1999).
131. Kurreck, J., Bieber, B., Jahn, R. & Erdmann, V. A. Comparative study of DNA enzymes and ribozymes against the same full-length messenger RNA of the vanilloid receptor subtype I. *J Biol Chem* **277**, 7099-107 (2002).
132. Grunweller, A. et al. Comparison of different antisense strategies in mammalian cells using locked nucleic acids, 2'-O-methyl RNA, phosphorothioates and small interfering RNA. *Nucleic Acids Res* **31**, 3185-93 (2003).
133. De Petrocellis, L. et al. The vanilloid receptor (VR1)-mediated effects of anandamide are potently enhanced by the cAMP-dependent protein kinase. *J Neurochem* **77**, 1660-3 (2001).

134. Gunthorpe, M. J., Harries, M. H., Prinjha, R. K., Davis, J. B. & Randall, A. Voltage- and time-dependent properties of the recombinant rat vanilloid receptor (rVR1). *J Physiol* **525 Pt 3**, 747-59 (2000).
135. Platika, D., Boulos, M. H., Baizer, L. & Fishman, M. C. Neuronal traits of clonal cell lines derived by fusion of dorsal root ganglia neurons with neuroblastoma cells. *Proc Natl Acad Sci U S A* **82**, 3499-503 (1985).
136. Francel, P. C. et al. Neurochemical characteristics of a novel dorsal root ganglion X neuroblastoma hybrid cell line, F-11. *J Neurochem* **48**, 1624-31 (1987).
137. Jahnel, R. et al. Dual expression of mouse and rat VRL-1 in the dorsal root ganglion derived cell line F-11 and biochemical analysis of VRL-1 after heterologous expression. *Eur J Biochem* **270**, 4264-71 (2003).
138. Olah, Z., Karai, L. & Iadarola, M. J. Anandamide activates vanilloid receptor 1 (VR1) at acidic pH in dorsal root ganglia neurons and cells ectopically expressing VR1. *J Biol Chem* **276**, 31163-70 (2001).
139. Bender, F. L. P. et al. Expression and functional characterization of the Vanilloid Receptor-like TRP channel VRL-1 in the primary sensory cell line F-11. *Abstract Viewer/Itinerary Planner. Washington, DC: Society for Neuroscience. Online.*, Program No. 48.23 (2002).
140. Goswami, C. et al. Identification and characterization of a Ca²⁺-sensitive interaction of the vanilloid receptor TRPV1 with tubulin. *J Neurochem* (2004).
141. Wilkins, M. R., Sanchez, J. C., Williams, K. L. & Hochstrasser, D. F. Current challenges and future applications for protein maps and post-translational vector maps in proteome projects. *Electrophoresis* **17**, 830-8 (1996).
142. Gygi, S. P., Rochon, Y., Franza, B. R. & Aebersold, R. Correlation between protein and mRNA abundance in yeast. *Mol Cell Biol* **19**, 1720-30 (1999).
143. Karas, M. & Hillenkamp, F. Laser desorption ionization of proteins with molecular masses exceeding 10,000 daltons. *Anal Chem* **60**, 2299-301 (1988).
144. Dresbach, T., Qualmann, B., Kessels, M. M., Garner, C. C. & Gundelfinger, E. D. The presynaptic cytomatrix of brain synapses. *Cell Mol Life Sci* **58**, 94-116 (2001).
145. Phillips, G. R. et al. The presynaptic particle web: ultrastructure, composition, dissolution, and reconstitution. *Neuron* **32**, 63-77 (2001).
146. Kennedy, M. B. The postsynaptic density. *Curr Opin Neurobiol* **3**, 732-7 (1993).
147. Gomperts, S. N. Clustering membrane proteins: It's all coming together with the PSD-95/SAP90 protein family. *Cell* **84**, 659-62 (1996).
148. Santoni, V., Molloy, M. & Rabilloud, T. Membrane proteins and proteomics: un amour impossible? *Electrophoresis* **21**, 1054-70 (2000).
149. Berggren, K. N. et al. An improved formulation of SYPRO Ruby protein gel stain: comparison with the original formulation and with a ruthenium II tris (bathophenanthroline disulfonate) formulation. *Proteomics* **2**, 486-98 (2002).
150. Gygi, S. P., Corthals, G. L., Zhang, Y., Rochon, Y. & Aebersold, R. Evaluation of two-dimensional gel electrophoresis-based proteome analysis technology. *Proc Natl Acad Sci U S A* **97**, 9390-5 (2000).
151. Godovac-Zimmermann, J. & Brown, L. R. Perspectives for mass spectrometry and functional proteomics. *Mass Spectrom Rev* **20**, 1-57 (2001).

152. Uetz, P. et al. A comprehensive analysis of protein-protein interactions in *Saccharomyces cerevisiae*. *Nature* **403**, 623-7 (2000).
153. Gavin, A. C. et al. Functional organization of the yeast proteome by systematic analysis of protein complexes. *Nature* **415**, 141-7 (2002).
154. Wang, H. & Hanash, S. Multi-dimensional liquid phase based separations in proteomics. *J Chromatogr B Analyt Technol Biomed Life Sci* **787**, 11-8 (2003).
155. Washburn, M. P., Wolters, D. & Yates, J. R., 3rd. Large-scale analysis of the yeast proteome by multidimensional protein identification technology. *Nat Biotechnol* **19**, 242-7 (2001).
156. Gygi, S. P. et al. Quantitative analysis of complex protein mixtures using isotope-coded affinity tags. *Nat Biotechnol* **17**, 994-9 (1999).
157. Zhou, H., Ranish, J. A., Watts, J. D. & Aebersold, R. Quantitative proteome analysis by solid-phase isotope tagging and mass spectrometry. *Nat Biotechnol* **20**, 512-5 (2002).
158. Amaya, F. et al. Local inflammation increases vanilloid receptor 1 expression within distinct subgroups of DRG neurons. *Brain Res* **963**, 190-6 (2003).
159. Tohda, C. et al. Axonal transport of VR1 capsaicin receptor mRNA in primary afferents and its participation in inflammation-induced increase in capsaicin sensitivity. *J Neurochem* **76**, 1628-35 (2001).
160. Carlton, S. M. & Coggeshall, R. E. Peripheral capsaicin receptors increase in the inflamed rat hindpaw: a possible mechanism for peripheral sensitization. *Neurosci Lett* **310**, 53-6 (2001).
161. Hudson, L. J. et al. VR1 protein expression increases in undamaged DRG neurons after partial nerve injury. *Eur J Neurosci* **13**, 2105-14 (2001).
162. Mevel-Ninio, M. & Weiss, M. C. Immunofluorescence analysis of the time-course of extinction, reexpression, and activation of albumin production in rat hepatoma-mouse fibroblast heterokaryons and hybrids. *J Cell Biol* **90**, 339-50 (1981).
163. Wainwright, A., Rutter, A. R., Seabrook, G. R., Reilly, K. & Oliver, K. R. Discrete expression of TRPV2 within the hypothalamo-neurohypophysial system: Implications for regulatory activity within the hypothalamic-pituitary-adrenal axis. *J Comp Neurol* **474**, 24-42 (2004).
164. Dietrich, A. et al. N-linked protein glycosylation is a major determinant for basal TRPC3 and TRPC6 channel activity. *J Biol Chem* **278**, 47842-52 (2003).
165. Watanabe, I. et al. Glycosylation affects rat Kv1.1 potassium channel gating by a combined surface potential and cooperative subunit interaction mechanism. *J Physiol* (2003).
166. Tyrrell, L., Renganathan, M., Dib-Hajj, S. D. & Waxman, S. G. Glycosylation alters steady-state inactivation of sodium channel Nav1.9/NaN in dorsal root ganglion neurons and is developmentally regulated. *J Neurosci* **21**, 9629-37 (2001).
167. Barnhill, J. C. et al. RGA protein associates with a TRPV ion channel during biosynthesis and trafficking. *J Cell Biochem* **91**, 808-20 (2004).
168. Hu, H. Z. et al. 2-Aminoethoxydiphenyl borate is a common activator of TRPV1, TRPV2, and TRPV3. *J Biol Chem* (2004).

169. Lazar, J., Szabo, T., Kovacs, L., Blumberg, P. M. & Biro, T. Distinct features of recombinant rat vanilloid receptor-1 expressed in various expression systems. *Cell Mol Life Sci* **60**, 2228-40 (2003).
170. Karai, L. J., Russell, J. T., Iadarola, M. J. & Olah, Z. Vanilloid receptor 1 regulates multiple calcium compartments and contributes to Ca²⁺-induced Ca²⁺ release in sensory neurons. *J Biol Chem* **279**, 16377-87 (2004).
171. Wisnoskey, B. J., Sinkins, W. G. & Schilling, W. P. Activation of vanilloid receptor type I in the endoplasmic reticulum fails to activate store-operated Ca²⁺ entry. *Biochem J* **372**, 517-28 (2003).
172. Guo, A., Vulchanova, L., Wang, J., Li, X. & Elde, R. Immunocytochemical localization of the vanilloid receptor 1 (VR1): relationship to neuropeptides, the P2X3 purinoceptor and IB4 binding sites. *Eur J Neurosci* **11**, 946-58 (1999).
173. Morenilla-Palao, C., Planells-Cases, R., Garcia-Sanz, N. & Ferrer-Montiel, A. Regulated exocytosis contributes to protein kinase C potentiation of vanilloid receptor activity. *J Biol Chem* **279**, 25665-72 (2004).
174. Bhawe, G. et al. Protein kinase C phosphorylation sensitizes but does not activate the capsaicin receptor transient receptor potential vanilloid 1 (TRPV1). *Proc Natl Acad Sci U S A* **100**, 12480-5 (2003).
175. Wang, C., Hu, H. Z., Colton, C. K., Wood, J. D. & Zhu, M. X. An alternative splicing product of the murine *trpv1* gene dominant-negatively modulates the activity of TRPV1 channels. *J Biol Chem* (2004).
176. Garcia-Sanz, N. et al. Identification of a tetramerization domain in the C terminus of the vanilloid receptor. *J Neurosci* **24**, 5307-14 (2004).
177. Watanabe, I., Zhu, J., Recio-Pinto, E. & Thornhill, W. B. Glycosylation affects the protein stability and cell surface expression of Kv1.4 but Not Kv1.1 potassium channels. A pore region determinant dictates the effect of glycosylation on trafficking. *J Biol Chem* **279**, 8879-85 (2004).
178. Keddi, N. et al. Analysis of the native quaternary structure of vanilloid receptor 1. *J Biol Chem* **276**, 28613-9 (2001).
179. Kuzhikandathil, E. V. et al. Functional analysis of capsaicin receptor (vanilloid receptor subtype 1) multimerization and agonist responsiveness using a dominant negative mutation. *J Neurosci* **21**, 8697-706 (2001).
180. Hofmann, T., Schaefer, M., Schultz, G. & Gudermann, T. Subunit composition of mammalian transient receptor potential channels in living cells. *Proc Natl Acad Sci U S A* **99**, 7461-6 (2002).
181. Pourrier, M., Schram, G. & Nattel, S. Properties, expression and potential roles of cardiac K⁺ channel accessory subunits: MinK, MiRPs, KChIP, and KChAP. *J Membr Biol* **194**, 141-52 (2003).
182. Ren, X., Shand, S. H. & Takimoto, K. Effective association of Kv channel-interacting proteins with Kv4 channel is mediated with their unique core peptide. *J Biol Chem* **278**, 43564-70 (2003).
183. Tyson, J. et al. IsK and KvLQT1: mutation in either of the two subunits of the slow component of the delayed rectifier potassium channel can cause Jervell and Lange-Nielsen syndrome. *Hum Mol Genet* **6**, 2179-85 (1997).
184. DePina, A. S. & Langford, G. M. Vesicle transport: the role of actin filaments and myosin motors. *Microsc Res Tech* **47**, 93-106 (1999).
185. Goldstein, L. S. & Yang, Z. Microtubule-based transport systems in neurons: the roles of kinesins and dyneins. *Annu Rev Neurosci* **23**, 39-71 (2000).

186. Guengerich, F. P. Common and uncommon cytochrome P450 reactions related to metabolism and chemical toxicity. *Chem Res Toxicol* **14**, 611-50 (2001).
187. Capdevila, J. H. & Falck, J. R. Biochemical and molecular properties of the cytochrome P450 arachidonic acid monooxygenases. *Prostaglandins Other Lipid Mediat* **68-69**, 325-44 (2002).
188. Kroetz, D. L. & Zeldin, D. C. Cytochrome P450 pathways of arachidonic acid metabolism. *Curr Opin Lipidol* **13**, 273-83 (2002).
189. Nedelkov, D. & Nelson, R. W. Delineation of in vivo assembled multiprotein complexes via biomolecular interaction analysis mass spectrometry. *Proteomics* **1**, 1441-6 (2001).
190. Velazquez-Campoy, A., Leavitt, S. A. & Freire, E. Characterization of protein-protein interactions by isothermal titration calorimetry. *Methods Mol Biol* **261**, 35-54 (2004).
191. Leavitt, S. & Freire, E. Direct measurement of protein binding energetics by isothermal titration calorimetry. *Curr Opin Struct Biol* **11**, 560-6 (2001).
192. Langton, P. E., Hankey, G. J. & Eikelboom, J. W. Cardiovascular safety of rofecoxib (Vioxx): lessons learned and unanswered questions. *Med J Aust* **181**, 524-5 (2004).
193. Sibbald, B. Rofecoxib (Vioxx) voluntarily withdrawn from market. *Cmaj* **171**, 1027-8 (2004).
194. Cortright, D. N. & Szallasi, A. Biochemical pharmacology of the vanilloid receptor TRPV1. An update. *Eur J Biochem* **271**, 1814-9 (2004).
195. Van Der Stelt, M. & Di Marzo, V. Endovanilloids. Putative endogenous ligands of transient receptor potential vanilloid 1 channels. *Eur J Biochem* **271**, 1827-34 (2004).
196. Heukeshoven, J. & Dernick, R. Increased sensitivity for Coomassie staining of sodium dodecyl sulfate-polyacrylamide gels using PhastSystem Development Unit. *Electrophoresis* **9**, 60-1 (1988).
197. Shevchenko, A., Wilm, M., Vorm, O. & Mann, M. Mass spectrometric sequencing of proteins silver-stained polyacrylamide gels. *Anal Chem* **68**, 850-8 (1996).
198. Saiki, R. K. et al. Enzymatic amplification of beta-globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* **230**, 1350-4 (1985).
199. Sambrook, J. & Russell, D. W. *Molecular Cloning: A laboratory manual* (Cold Spring Laboratory Press, Cold Spring Harbor, New York, 2001).
200. Hanahan, D. Studies on transformation of *Escherichia coli* with plasmids. *J Mol Biol* **166**, 557-80 (1983).
201. Birnboim, H. C. & Doly, J. A rapid alkaline extraction procedure for screening recombinant plasmid DNA. *Nucleic Acids Res* **24**, 1513-23 (1979).
202. Sanger, F., Nicklen, S. & Coulson, A. R. DNA sequencing with chain-terminating inhibitors. *Proc Natl Acad Sci U S A* **74**, 5463-7 (1977).
203. Laemmli, U. K. Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature* **227**, 680-5 (1970).
204. Schagger, H. & von Jagow, G. Blue native electrophoresis for isolation of membrane protein complexes in enzymatically active form. *Anal Biochem* **199**, 223-31 (1991).

205. Schagger, H., Cramer, W. A. & von Jagow, G. Analysis of molecular masses and oligomeric states of protein complexes by blue native electrophoresis and isolation of membrane protein complexes by two-dimensional native electrophoresis. *Anal Biochem* **217**, 220-30 (1994).
206. Nadano, D., Aoki, C., Yoshinaka, T., Irie, S. & Sato, T. A. Electrophoretic characterization of ribosomal subunits and proteins in apoptosis: specific downregulation of S11 in staurosporine-treated human breast carcinoma cells. *Biochemistry* **40**, 15184-93 (2001).
207. Brown, R. S. & Lennon, J. J. Mass resolution improvement by incorporation of pulsed ion extraction in a matrix-assisted laser desorption/ionization linear time-of-flight mass spectrometer. *Anal Chem* **67**, 1998-2003 (1995).
208. Whittall, R. M. & Li, L. High-resolution matrix-assisted laser desorption/ionization in a linear time-of-flight mass spectrometer. *Anal Chem* **67**, 1950-4 (1995).
209. Spengler, B., Kirsch, D., Kaufmann, R. & Jaeger, E. Peptide sequencing by matrix-assisted laser-desorption mass spectrometry. *Rapid Commun Mass Spectrom* **6**, 105-8 (1992).
210. Kaufmann, R., Spengler, B. & Lutzenkirchen, F. Mass spectrometric sequencing of linear peptides by product-ion analysis in a reflectron time-of-flight mass spectrometer using matrix-assisted laser desorption ionization. *Rapid Commun Mass Spectrom* **7**, 902-10 (1993).